
```

% Copyright © 2019 Naturalpoint
%   Not for Redistribution
%   NatNet Matlab Class
%   Requires Matlab 2019a or later
%   This class is a wrapper for the NatNetML assembly and controls the
%   connection and communication between the OptiTrack NatNet Server and
%   Matlab clients.
%   Public Properties:
%       ClientIP - The IP address of the client network.
%       ConnectionType - Multicast or Unicast
%       HostIP - The IP address of the host/server network.
%       IsReporting - state of commandline reporting
%       IsConnected - state of connection to server, updated only when
%                   establishing a connection. (Read-only)
%       Mode - state of the Motive server, Live or Edit. (Read-Only)
%       Version - version number of the natnet libraries. (Read-Only)
%       FrameRate - sampling rate of the Motive hardware. (Read-Only)
%   Public Methods:
%       natnet - Object Constructor
%       connect - Establish a connection to an OptiTrack NatNet Server
%       disconnect - Close a connection to an OptiTrack NatNet Server
%       addlistener - Links a callback function to the server stream.

%       getFrameRate - Returns the server Frame Rate.
%       getModelDescription - Returns a structure of the asset list.
%       getFrame - Returns the latest frame of mocap data.

%       getFrameMetaData - Returns last frame of meta data within mocap
%                           frame. Return value is a Matlab struct
%       getFrameMarkerSet - Returns the last frame of MarkerSet data
%                           within mocap frame. Return value is a Matlab struct
%       getFrameLabeledMarker - Returns the last frame of LabeledMarker
%                               data within the mocap frame. Return value is a Matlab
%                               struct.

%       getMode - Returns the state of the server: Live or Edit mode.
%       enable - Enables execution of a linked callback function.
%       disable - Disabled execution of a linked callback function.
%       Server Commands
%       startRecord - Starts a recording on the server.
%       stopRecord - Stops a recording on the server.
%       cycleRecord - Loops recording.
%       liveMode - Changes the server to Live mode.
%       editMode - Changes the server to Edit mode.
%       startPlayback - Starts playback of the server in Edit mode.
%       stopPlayback - Stops playback of the server in Edit mode.
%       setPlaybackStartFrame - Sets the playback loop start frame.
%       setPlaybackEndFrame - Sets the playback loop end frame.
%       setPlaybackLooping - Sets playback to loop.
%       setPlaybackCurrentFrame - Sets the current frame in Edit mode.
%       setTakeName - Sets the pending (to be recorded) take name.
%       setPlaybackTakeName - Opens a recorded take in Edit mode.

```

```

%      setCurrentSession - Changes the current folder in Edit mode.
%      delete - Destructor for the natnet class

classdef natnet < handle
    properties ( Access = public )
        ClientIP
        ConnectionType
        HostIP
        IsReporting
    end % properties ( Access = public )

    properties ( SetAccess = private )
        IsConnected
        Mode
        Version
        FrameRate
    end % properties ( SetAccess = private )

    properties ( Access = private )
        Assembly
        AssemblyPath
        AssemblyVersion
        Client
        IsAssemblyLoaded
        LastAssemblyPathFile
        MaxListeners
        Listener
        CReattempt
    end % properties ( Access = private )

    properties ( Access = private , Dependent = true )
        iConnectionType
    end % properties ( Access = private , Dependant = true )

    methods
        function set.HostIP( obj , val )
            validIP = obj.checkip( val);
            if ( ~isempty( validIP ) )
                obj.disconnect
                obj.HostIP = validIP;
                obj.report( [ 'set host ip address: ' , validIP ] )
            end
        end % set.HostIP

        function set.ClientIP( obj , val )
            validIP = obj.checkip( val );
            if ( ~isempty( validIP ) )
                obj.disconnect
                obj.ClientIP = validIP;
            end
        end % set.ClientIP
    end
end

```

```

        obj.report( [ 'set client ip address: ' , validIP ] )
    end
end % set.ClientIP

function set.ConnectionType( obj , val )
    obj.disconnect
    if ( strcmpi( 'Multicast' , val ) == 1 )
        obj.ConnectionType = 'Multicast';
        obj.report( 'set multicast' )
    elseif ( strcmpi( 'Unicast' , val ) == 1 )
        obj.ConnectionType = 'Unicast';
        obj.report( 'set unicast' )
    else
        obj.report( 'invalid connection type' )
    end
end % set.ConnectionType

function set.IsReporting( obj , val )
    if ( val == 1 )
        obj.IsReporting = 1;
        obj.report( 'reporting enabled' )
    elseif ( val == 0 )
        obj.report( 'reporting disabled' )
        obj.IsReporting = 0;
    end
end % set.IsReporting

function iconnectiontype = get.iConnectionType( obj )
    if ( strcmpi( 'Multicast' , obj.ConnectionType ) == 1 )
        iconnectiontype = 0;
    elseif ( strcmpi( 'Unicast' , obj.ConnectionType ) == 1 )
        iconnectiontype = 1;
    end
end % iconnectiontype
end % methods ( property set and get methods)

methods ( Access = public )
    function obj = natnet( )
        obj.getLastAssemblyPath
        obj.MaxListeners = 255;
        obj.Listener{ obj.MaxListeners , 1 } = [ ];
        obj.FrameRate = 0;
        obj.IsAssemblyLoaded = 0;
        obj.IsConnected = 0;
        obj.IsReporting = 0;
        obj.CReattempt = 1;
        obj.HostIP = '127.0.0.1';
        obj.ClientIP = '127.0.0.1';
        obj.ConnectionType = 'Multicast';
        obj.Version = ' ';
    end
end

```

```

        obj.Mode = ' ';
end % natnet - constructor

function connect( obj )
    obj.disconnect
    obj.getAssemblies
    if ( obj.IsAssemblyLoaded == 0 )
        if ( isempty( obj.AssemblyPath ) )
            obj.setAssemblyPath
            obj.addAssembly
            obj.getAssemblies
        else
            obj.addAssembly
            obj.getAssemblies
        end
    end
    if ( obj.IsAssemblyLoaded == 0 )
        report( obj , 'natnetml assembly is missing or undefined' )
        return
    end
    try
        obj.Client = NatNetML.NatNetClientML( obj.iConnectionType );
        v = obj.Client.NatNetVersion( );
        obj.Version = sprintf( '%d.%d.%d.%d', v( 1 ) , v( 2 ) ,
v( 3 ) , v( 4 ) );
        if ( isempty( obj.Client ) == 1 || obj.IsAssemblyLoaded ==
0 )
            obj.report( 'client object invalid' )
            return
        else
            obj.report( [ 'natnet version: ' , obj.Version ] )
        end
        flg = obj.Client.Initialize( obj.ClientIP , obj.HostIP );
        if flg == 0
            obj.report( 'client initialized' )
        else
            obj.report( 'initialization failed' )
            obj.IsConnected = 0;
            return
        end
    catch exception
        rethrow(exception)
    end

    [ ~ , rc] = obj.Client.SendMessageAndWait( 'FrameRate' );
    if rc == 0
        obj.IsConnected = 1;
        obj.getFrameRate
        obj.getMode
        obj.getModelDescription;
        return
    elseif( rc == 1 )
        obj.report( 'connection failed due to an internal error' )

```

```

elseif ( rc == 2 )
    obj.report( 'connection failed due to an external error' )
elseif ( rc == 3 )
    obj.report( 'connection failed due to a network error' )
else
    obj.report( 'connection failed due to an unknown error' )
end
obj.IsConnected = 0;
end % connect

function addlistener( obj , aListenerIndex , functionhandlestring )
    if ~isnumeric( aListenerIndex ) || aListenerIndex < 1 ||
aListenerIndex > obj.MaxListeners
        obj.report( 'invalid index' );
        return
    end

    if isa( functionhandlestring , 'char' ) &&
~isempty( which( functionhandlestring ) )
        functionhandle = str2func ( [ '@( src , evnt)' ,
functionhandlestring , '( src , evnt)' ] );

        if ~isa( obj.Listener{ aListenerIndex , 1 } ,
'event.listener' )
            obj.report( [ 'adding listener in slot:' ,
num2str(aListenerIndex ) ] );
            obj.Listener{ aListenerIndex , 1 } =
addlistener( obj.Client , 'OnFrameReady2' , functionhandle );
            obj.disable( aListenerIndex );
        elseif isa( obj.Listener{ aListenerIndex , 1 } ,
'event.listener' )
            obj.report( [ 'replacing listener in slot: ' ,
num2str(aListenerIndex ) ] )
            obj.disable( aListenerIndex )
            delete( obj.Listener{ aListenerIndex , 1 } )
            obj.Listener{ aListenerIndex , 1 } =
addlistener( obj.Client , 'OnFrameReady2' , functionhandle );
            obj.disable ( aListenerIndex )
        end
    else
        obj.report( [ 'invalid function set for slot: ' ,
num2str( aListenerIndex ) ] )
        return
    end
end % addlistener

function getFrameRate( obj )
    [ bytearray , rc ] = obj.sendMessageAndWait( 'FrameRate' );
    if rc == 0
        bytearray = uint8( bytearray );
        obj.FrameRate = typecast( bytearray , 'single' );
        obj.report( [ 'frame rate: ' , num2str( obj.FrameRate ) , '

```

```

    fps' ] )
        else
            obj.FrameRate = 0;
        end
    end % getFrameRate

    function modelDescription = getModelDescription( obj )
        if ( obj.IsConnected == 1 )
            dataDescriptions = obj.Client.GetDataDescriptions( );
            modelDescription.TrackingModelCount = dataDescriptions.Count;
            report( obj , [ 'number of trackables: ' ,
num2str( dataDescriptions.Count ) ] )
            modelDescription.MarkerSetCount = 0;
            modelDescription.RigidBodyCount = 0;
            modelDescription.SkeletonCount = 0;
            modelDescription.ForcePlateCount = 0;
            modelDescription.DeviceCount = 0;
            modelDescription.CameraCount = 0;

            for i = 1 : modelDescription.TrackingModelCount
                descriptor = dataDescriptions.Item( i - 1 );
                % marker sets
                if ( descriptor.type == 0 )
                    modelDescription.MarkerSetCount =
modelDescription.MarkerSetCount + 1;

modelDescription.MarkerSet( modelDescription.MarkerSetCount ).Name =
char( descriptor.Name );

modelDescription.MarkerSet( modelDescription.MarkerSetCount ).MarkerCount =
descriptor.nMarkers;
                    for k = 1 : descriptor.nMarkers

modelDescription.MarkerSet( modelDescription.MarkerSetCount ).Markers( k ).La
bel = char( descriptor.MarkerNames( k ) );
                        end
                        % rigid bodies
                        elseif ( descriptor.type == 1 )
                            modelDescription.RigidBodyCount =
modelDescription.RigidBodyCount + 1;

modelDescription.RigidBody( modelDescription.RigidBodyCount ).Name =
char( descriptor.Name );

modelDescription.RigidBody( modelDescription.RigidBodyCount ).ID =
descriptor.ID ;

modelDescription.RigidBody( modelDescription.RigidBodyCount ).ParentID =
descriptor.parentID;

modelDescription.RigidBody( modelDescription.RigidBodyCount ).OffsetX =
descriptor.offsetx;

```

```

modelDescription.RigidBody( modelDescription.RigidBodyCount ).OffsetX =
descriptor.offsetx;

modelDescription.RigidBody( modelDescription.RigidBodyCount ).OffsetX =
descriptor.offsetx;

modelDescription.RigidBody( modelDescription.RigidBodyCount ).Type =
descriptor.type;
    % skeletons
    elseif ( descriptor.type == 2 )
        modelDescription.SkeletonCount =
modelDescription.SkeletonCount + 1;

modelDescription.Skeleton( modelDescription.SkeletonCount ).Name =
char( descriptor.Name );

modelDescription.Skeleton( modelDescription.SkeletonCount ).ID =
descriptor.ID;

modelDescription.Skeleton( modelDescription.SkeletonCount ).SegmentCount =
descriptor.nRigidBodies;
        for k = 1 :
modelDescription.Skeleton( modelDescription.SkeletonCount ).SegmentCount

modelDescription.Skeleton( modelDescription.SkeletonCount ).Segment( k ).Name
    = char( descriptor.RigidBodies( k ).Name );

modelDescription.Skeleton( modelDescription.SkeletonCount ).Segment( k ).ID
= descriptor.RigidBodies( k ).ID;

modelDescription.Skeleton( modelDescription.SkeletonCount ).Segment( k ).ParentID = descriptor.RigidBodies( k ).parentID;

modelDescription.Skeleton( modelDescription.SkeletonCount ).Segment( k ).OffsetX = descriptor.RigidBodies( k ).offsetx;

modelDescription.Skeleton( modelDescription.SkeletonCount ).Segment( k ).OffsetY = descriptor.RigidBodies( k ).offsety;

modelDescription.Skeleton( modelDescription.SkeletonCount ).Segment( k ).OffsetZ = descriptor.RigidBodies( k ).offsetz;

modelDescription.Skeleton( modelDescription.SkeletonCount ).Segment( k ).Type
    = descriptor.RigidBodies( k ).type;
        end
        elseif ( descriptor.type == 3 )
            modelDescription.ForcePlateCount =
modelDescription.ForcePlateCount + 1;

modelDescription.ForcePlate( modelDescription.ForcePlateCount ).Serial =
char( descriptor.Serial );

modelDescription.ForcePlate( modelDescription.ForcePlateCount ).ID =
descriptor.ID;

```

```

modelDescription.ForcePlate( modelDescription.ForcePlateCount ).Width =
descriptor.Width;

modelDescription.ForcePlate( modelDescription.ForcePlateCount ).Length =
descriptor.Length;

modelDescription.ForcePlate( modelDescription.ForcePlateCount ).Origin.X =
descriptor.OriginX;

modelDescription.ForcePlate( modelDescription.ForcePlateCount ).Origin.Y =
descriptor.OriginY;

modelDescription.ForcePlate( modelDescription.ForcePlateCount ).Origin.Z =
descriptor.OriginZ;

modelDescription.ForcePlate( modelDescription.ForcePlateCount ).CalibrationMa
trix = descriptor.CalMatrix;

modelDescription.ForcePlate( modelDescription.ForcePlateCount ).CornerCount
= descriptor.Corners.Length;
    for i = 1 :
modelDescription.ForcePlate( modelDescription.ForcePlateCount ).CornerCount/3

modelDescription.ForcePlate( modelDescription.ForcePlateCount ).Corner( i ).X
= descriptor.Corners( i*3-2 );

modelDescription.ForcePlate( modelDescription.ForcePlateCount ).Corner( i ).Y
= descriptor.Corners( i*3-1 );

modelDescription.ForcePlate( modelDescription.ForcePlateCount ).Corner( i ).Z
= descriptor.Corners( i*3 );
    end

modelDescription.ForcePlate( modelDescription.ForcePlateCount ).PlateType =
descriptor.PlateType;

modelDescription.ForcePlate( modelDescription.ForcePlateCount ).ChannelDataTy
pe = descriptor.ChannelDataType;

modelDescription.ForcePlate( modelDescription.ForcePlateCount ).ChannelCount
= descriptor.ChannelCount;
    for i = 1 :
modelDescription.ForcePlate( modelDescription.ForcePlateCount ).ChannelCount

modelDescription.ForcePlate( modelDescription.ForcePlateCount ).Channel( i ).
Name = char( descriptor.ChannelNames( i ) );
    end
    elseif (descriptor.type == 4)
        modelDescription.DeviceCount =
modelDescription.DeviceCount + 1;
    elseif (descriptor.type == -1 || descriptor.type == 5)
        modelDescription.CameraCount =
modelDescription.CameraCount + 1;

```

```

        else
            report( obj , 'invalid asset type' )
        end
    end
    report( obj , [ 'number of markersets: ' ,
num2str( modelDescription.MarkerSetCount ) ] )
    report( obj , [ 'number of rigid bodies: ' ,
num2str( modelDescription.RigidBodyCount ) ] )
    report( obj , [ 'number of skeletons: ' ,
num2str( modelDescription.SkeletonCount ) ] )
    report( obj , [ 'number of force plates: ' ,
num2str( modelDescription.ForcePlateCount ) ] )
    report( obj , [ 'number of devices: ' ,
num2str( modelDescription.DeviceCount ) ] )
    report( obj , [ 'number of cameras: ' ,
num2str( modelDescription.CameraCount ) ] )
    else
        obj.report( 'connection not established' )
    end
end % getModelDescription

function frameOfMocapData = getFrame( obj )
    if ( obj.IsConnected == 1 )
        data = obj.Client.GetLastFrameOfData( );
        frameOfMocapData = data;
    else
        obj.report( 'connection not established' )
    end
end % getFrame

function frameOfMetaData = getFrameMetaData( obj )
    persistent data
    if ( obj.IsConnected == 1 )
        data = obj.Client.GetLastFrameOfData( );
        t.Frame = data.iFrame;
        t.Timestamp = data.fTimestamp;
        t.CameraMidExposureTimestamp =
data.CameraMidExposureTimestamp;
        t.CameraDataReceivedTimestamp =
data.CameraDataReceivedTimestamp;
        t.TransmitTimestamp = data.TransmitTimestamp;
        t.Recording = data.bRecording;
        t.TrackingModelsChanged = data.bTrackingModelsChanged;
        t.Timecode = data.Timecode;
        t.TimecodeSubframe = data.TimecodeSubframe;
        t.MarkerCount = data.nMarkers;
        t.OtherMarkerCount = data.nOtherMarkers;
        t.MarkerSetCount = data.nMarkerSets;
        t.RigidBodyCount = data.nRigidBodies;
        t.SkeletonCount = data.nSkeletons;
        t.ForcePlateCount = data.nForcePlates;
        t.DeviceCount = data.nDevices;
    end
end % getFrameMetaData

```

```

        frameOfMetaData = t;
    else
        obj.report( 'connection not established' )
    end
end % getFrameOfMetaData

function frameOfLabeledMarkerData = getFrameLabeledMarker( obj )
    if ( obj.IsConnected == 1 )
        data = obj.Client.GetLastFrameOfData( );
        t.Frame = data.iFrame;
        t.Timestamp = data.fTimestamp;
        t.MarkerCount = data.nMarkers;
        for i = t.MarkerCount : -1 : 1
            t.LabeledMarker( i ).AssetID =
                int16( bitshift( data.LabeledMarkers( i ).ID , -16 ) );
            t.LabeledMarker( i ).MemberID = int16( bitshift
                ( bitshift( data.LabeledMarkers( i ).ID , 16 ) , -16 ) );
            t.LabeledMarker( i ).X = data.LabeledMarkers( i ).x;
            t.LabeledMarker( i ).Y = data.LabeledMarkers( i ).y;
            t.LabeledMarker( i ).Z = data.LabeledMarkers( i ).z;
            t.LabeledMarker( i ).Size =
                data.LabeledMarkers( i ).size;
            %t.LabeledMarker( i ).Occluded =
                logical( bitget( data.LabeledMarkers( i ).parameters , 1 ) );
            t.LabeledMarker( i ).PointCloudSolved =
                logical( bitget( data.LabeledMarkers( i ).parameters , 2 ) );
            t.LabeledMarker( i ).AssetSolved =
                logical( bitget( data.LabeledMarkers( i ).parameters , 3 ) );
            t.LabeledMarker( i ).AssetMarker =
                logical( bitget( data.LabeledMarkers( i ).parameters , 4 ) );
            t.LabeledMarker( i ).UnlabeledMarker =
                logical( bitget( data.LabeledMarkers( i ).parameters , 5 ) );
            t.LabeledMarker( i ).ActiveMarker =
                logical( bitget( data.LabeledMarkers( i ).parameters , 6 ) );
            t.LabeledMarker( i ).Residual =
                data.LabeledMarkers( i ).residual;
        end
        frameOfLabeledMarkerData = t;
    else
        obj.report( 'connection not established' )
    end
end % getFrameLabeledMarker

function frameOfMarkerSetData = getFrameMarkerSet( obj )
    if ( obj.IsConnected == 1 )
        data = obj.Client.GetLastFrameOfData( );
        t.Frame = data.iFrame;
        t.Timestamp = data.fTimestamp;
        t.MarkerSetCount = data.nMarkerSets;
        for i = t.MarkerSetCount : -1 : 1
            t.MarkerSet( i ).MarkerSetName =
                string( data.MarkerSets( i ).MarkerSetName );
        end
    end
end % getFrameMarkerSet

```

```

        t.MarkerSet( i ).MarkerCount =
data.MarkerSets( i ).nMarkers;
        for j = t.MarkerSet( i ).MarkerCount : -1 : 1
            %t.MarkerSet( i ).Marker( j ).AssetID =
intl6( bitshift( data.MarkerSets( i ).Markers( j ).ID , -16 ) );
            %t.MarkerSet( i ).Marker( j ).MemberID =
intl6( bitshift ( bitshift( data.MarkerSets( i ).Markers( j ).ID , 16 ) ,
-16 ) );
            t.MarkerSet( i ).Marker( j ).X =
data.MarkerSets( i ).Markers( j ).x;
            t.MarkerSet( i ).Marker( j ).Y =
data.MarkerSets( i ).Markers( j ).y;
            t.MarkerSet( i ).Marker( j ).Z =
data.MarkerSets( i ).Markers( j ).z;
            %t.MarkerSet( i ).Marker( j ).Size =
data.MarkerSets( i ).Markers( j ).size;
            %t.MarkerSet( i ).Marker( j ).Occluded =
logical( bitget( data.MarkerSets( i ).Markers( j ).parameters , 1 ) );
            %t.MarkerSet( i ).Marker( j ).PointCloudSolved =
logical( bitget( data.MarkerSets( i ).Markers( j ).parameters , 2 ) );
            %t.MarkerSet( i ).Marker( j ).AssetSolved =
logical( bitget( data.MarkerSets( i ).Markers( j ).parameters , 3 ) );
            %t.MarkerSet( i ).Marker( j ).AssetMarker =
logical( bitget( data.MarkerSets( i ).Markers( j ).parameters , 4 ) );
            %t.MarkerSet( i ).Marker( j ).UnlabeledMarker =
logical( bitget( data.MarkerSets( i ).Markers( j ).parameters , 5 ) );
            %t.MarkerSet( i ).Marker( j ).ActiveMarker =
logical( bitget( data.MarkerSets( i ).Markers( j ).parameters , 6 ) );
            %t.MarkerSet( i ).Marker( j ).Residual =
data.MarkerSets( i ).Markers( j ).residual;
        end
    end
    frameOfMarkerSetData = t;
else
    obj.report( 'connection not established' )
end
end % setFrameMarkerSet

function enable( obj , eListenerIndex )
    if ~isnumeric( eListenerIndex ) || eListenerIndex < 0 ||
eListenerIndex > obj.MaxListeners
        obj.report( 'invalid index' )
        return
    end

    if obj.IsConnected == 0
        return
    end

    if eListenerIndex == 0
        for k = 1:obj.MaxListeners
            if isa( obj.Listener{ k , 1} , 'event.listener' ) &&
obj.Listener{ k , 1}.Enabled == false

```

```

        obj.Listener{ k , 1 }.Enabled = true;
        obj.report( [ 'listener enabled in slot: ' ,
num2str( k ) ] )
    end
end
else
    if( isa( obj.Listener{ eListenerIndex , 1 } ,
'event.listener' ) ) && obj.Listener{ eListenerIndex , 1 }.Enabled == false
        obj.Listener{ eListenerIndex , 1 }.Enabled = true;
        obj.report( [ 'listener enabled in slot: ' ,
num2str( eListenerIndex ) ] )
    end
end
end % enable

function disable( obj , dListenerIndex );
    if ~isnumeric( dListenerIndex ) || dListenerIndex < 0 ||
dListenerIndex > obj.MaxListeners
        obj.report( 'invalid index' );
        return
    end

    if obj.IsConnected == 0;
        return
    end

    if dListenerIndex == 0
        for k = 1:obj.MaxListeners
            if isa( obj.Listener{ k , 1 } , 'event.listener' ) &&
obj.Listener{ k , 1 }.Enabled == true;
                obj.Listener{ k }.Enabled = false;
                obj.report( [ 'listener disabled in slot: ' ,
num2str( k ) ] )
            end
        end
    else
        if( isa( obj.Listener{ dListenerIndex , 1 } ,
'event.listener' ) ) && obj.Listener{ dListenerIndex , 1 }.Enabled == true;
            obj.Listener{ dListenerIndex , 1 }.Enabled = false;
            obj.report( [ 'listener disabled in slot: ' ,
num2str( dListenerIndex ) ] )
        end
    end
end % disable

function startRecord( obj )
    obj.sendMessageAndWait( 'StartRecording' );
end % startRecord

function stopRecord( obj )
    obj.sendMessageAndWait( 'StopRecording' );

```

```

    end % stopRecord

    function cycleRecord( obj , iterations , duration , delay )
        if ( isnumeric( iterations ) && isnumeric( duration ) &&
isnumeric( delay ) );
            for i = 1:iterations
                pause( delay );
                obj.startRecord;
                pause( duration );
                obj.stopRecord;
            end
        end
    end % cycleRecord

    function liveMode( obj )
        obj.sendMessageAndWait( 'LiveMode' );
    end % liveMode

    function editMode( obj )
        obj.sendMessageAndWait( 'EditMode' );
    end % editMode

    function startPlayback( obj )
        obj.sendMessageAndWait( 'TimelinePlay' );
    end % startPlayback

    function stopPlayback( obj )
        obj.sendMessageAndWait( 'TimelineStop' );
    end % stopPlayback

    function setPlaybackStartFrame( obj , startFrame )
        obj.sendMessageAndWait( [ 'SetPlaybackStartFrame,' ,
num2str( startFrame ) ] );
    end % setPlaybackStartFrame

    function setPlaybackEndFrame( obj , endFrame )
        obj.sendMessageAndWait( [ 'SetPlaybackStopFrame,' ,
num2str( endFrame ) ] );
    end % setPlaybackEndFrame

    function setPlaybackLooping( obj , val )
        obj.sendMessageAndWait( [ 'SetPlaybackLooping,' ,
num2str( val ) ] );
    end % stopPlaybackLooping

```

```

function setPlaybackCurrentFrame( obj , currentFrame )
    obj.sendMessageAndWait( [ 'SetPlaybackCurrentFrame,' ,
num2str( currentFrame ) ] );
end % stopPlaybackLooping

function setTakeName( obj , name )
    obj.sendMessageAndWait( strcat( 'SetRecordTakeName,' , name ) );
end % setTakeName

function setPlaybackTakeName( obj , name )
    obj.sendMessageAndWait( strcat( 'SetPlaybackTakeName,' ,
name ) );
end % setPlaybackTakeName

function setCurrentSession( obj , name )
    obj.sendMessageAndWait( strcat( 'SetCurrentSession,' , name ) );
end % setCurrentSession

function getMode( obj )
    [ bytearray , rc ] = obj.sendMessageAndWait( 'CurrentMode' );
    if rc == 0
        state = bytearray(1);
        if ( state == 0 )
            obj.Mode = 'Live';
        elseif ( state == 1 )
            obj.Mode = 'Recording';
        elseif ( state == 2 )
            obj.Mode = 'Edit';
        else
            obj.Mode = ' ';
        end
        obj.report( [ 'mode: ' , lower( obj.Mode ) ] )
        return
    else
        obj.Mode = ' ';
    end
end % getServerState

function disconnect( obj )
    obj.disable( 0 );
    for k = 1 : obj.MaxListeners
        if isa( obj.Listener{ k , 1 } , 'event.listener' )
            delete( obj.Listener{ k , 1 } )
        end
    end
    obj.Listener{ obj.MaxListeners , 1 } = [ ];
    if ( isempty( obj.Client ) == 0 )
        obj.Client.Uninitialize;
        obj.Client = [ ];
    end
end

```

```

        if ( obj.IsConnected == 1 )
            report( obj , 'disconnecting' );
        end
        obj.IsConnected = 0;
    end
end % disconnect

function delete( obj )
    obj.disconnect
end % delete - destructor
end % methods ( Access = public )

methods ( Access = private )
    function getLastAssemblyPath( obj )
        pathtomfile = mfilename( 'fullpath' );
        mfile = mfilename( );
        mfilenamelength = length( mfile );
        foldertomfile = pathtomfile( 1 : end-mfilenamelength );
        obj.LastAssemblyPathFile = strcat( foldertomfile ,
'assemblypath.txt' );
        if ( exist( obj.LastAssemblyPathFile , 'file' ) == 2 )
            assemblypath = textread( obj.LastAssemblyPathFile , '%s' );
            assemblypath = strjoin( assemblypath );
            if( exist( assemblypath , 'file' ) == 2)
                obj.AssemblyPath = assemblypath;
            end
        end
    end
end % getLastAssemblyPath

function getAssemblies( obj )
    obj.IsAssemblyLoaded = 0;
    obj.Assembly = [ ];
    obj.AssemblyVersion = [ ];
    domain = System.AppDomain.CurrentDomain;
    assemblies = domain.GetAssemblies;
    assembly{ assemblies.Length , 1 } = [ ];
    for i = 1:assemblies.Length
        assembly{ i } = assemblies.Get( i-1 );
        obj.Assembly{ i } = char( assembly{ i }.FullName );
        if ( strncmp( obj.Assembly{ i } , 'NatNetML' , 8 ) )
            aver = regexp( obj.Assembly{ i } , '\d+' , 'match' );
            obj.AssemblyVersion = strcat( aver{ 1 } , '.' ,
aver{ 2 } , '.' , aver{ 3 } , '.' , aver{ 4 } );
            obj.IsAssemblyLoaded = 1;
        end
    end
end % getAssemblies

function addAssembly( obj )
    if ( exist( obj.AssemblyPath , 'file' ) == 2 )

```

```

        [ ~ ] = NET.addAssembly( obj.AssemblyPath );
        obj.IsAssemblyLoaded = 1;
    else
        obj.IsAssemblyLoaded = 0;
        report( obj , 'failed to add NatNetML.dll assembly' );
    end
end % addAssembly

function setAssemblyPath( obj )
    [ name , path ] = uigetfile( '*.dll' , 'Select the
NatNetML.dll assembly - NatNetLib.dll is a dependency' );
    if ( strcmpi( name , 'NatNetML.dll' ) == 1 )
        assemblyPath = strcat( path , name );
        textAssemblyPath = strrep( assemblyPath , '\\' , '\\\' );
        fileid = fopen( obj.LastAssemblyPathFile , 'wt' );
        fprintf( fileid , '%s' , assemblyPath );
        fclose( fileid );
        if ( obj.IsAssemblyLoaded == 0 )
            obj.AssemblyPath = assemblyPath;
            obj.report( [ 'defined assembly path: ' ,
textAssemblyPath ] )
        else
            obj.report( [ 'redefined assembly path: ' ,
textAssemblyPath ] )
            obj.report( 'restart matlab to apply the changes' )
        end
    else
        obj.report( 'invalid assembly path in function
setAssemblyPath' )
    end
end % setAssemlbyPath

function report( obj , message )
    if ( obj.IsReporting == 1 )
        ctime = strsplit( num2str( clock ) );
        disp( sprintf ( [ ctime{ 1 } , '/' , ctime{ 2 } , '/' ,
ctime{ 3 } , ' ' , ctime{ 4 } , ':' , ctime{ 5 } , ':' , ctime{ 6 } , '
' , 'natnet - ' , message ] ) )
    end
end % report

function validIP = checkip( obj , value )
    if ( ischar( value ) && length( value ) < 16 && length( value )
> 6 )
        val = strsplit( value , '.' );
        if ( length( val ) == 4 )
            if all( isstrprop( val{ 1 } , 'digit' ) ) &&
all( isstrprop( val{ 2 } , 'digit' ) )...
                && all( isstrprop( val{ 3 } , 'digit' ) ) &&
all( isstrprop( val{ 4 } , 'digit' ) )...
                && length( val{ 1 } ) < 4 && length( val{ 2 } )

```

```

< 4 && length( val{ 3 } ) < 4 && length( val{ 4 } ) < 4 ...
                                && ~isempty( val{ 1 } ) && ~isempty( val{ 2 } )
&& ~isempty( val{ 3 } ) && ~isempty( val{ 4 } ) ...
                                && str2double( val{ 1 } ) < 256 &&
str2double( val{ 2 } ) < 256 && str2double( val{ 3 } ) < 256 &&
str2double( val{ 4 } ) < 256 ...
                                && str2double( val{ 1 } ) >= 0 &&
str2double( val{ 2 } ) >= 0 && str2double( val{ 3 } ) >= 0 &&
str2double( val{ 4 } ) >= 0
                                validIP = value;
                                else
                                    report( obj , 'invalid string for ip address (e.x.
127.0.0.1)' )
                                    validIP = [ ];
                                end
                                else
                                    report( obj , 'invalid string for ip address (e.x.
127.0.0.1)' )
                                    validIP = [ ];
                                end
                                else
                                    report( obj , 'invalid string for ip address (e.x.
127.0.0.1)' )
                                    validIP = [ ];
                                end
                                end % validIP

function [ bytearray , rc ] = sendMessageAndWait( obj , cmd )
    if ( obj.IsConnected == 1 );
        [ bytearray , rc ] = obj.Client.SendMessageAndWait( cmd );
        if ( rc == 0 )
            return
        else
            obj.report( 'command failed due to an unknown error' )
        end
    else
        bytearray = '';
        rc = '';
        obj.report( 'connection not established' )
    end
end % sendMessageAndWait
end % methods ( Access = private )
end % classdef natnet < handle

ans =

    natnet with properties:

        ClientIP: '127.0.0.1'
        ConnectionType: 'Multicast'
        HostIP: '127.0.0.1'
        IsReporting: 0

```

IsConnected: 0
Mode: ' '
Version: ' '
FrameRate: 0

Published with MATLAB® R2023b